

Catalog Number

9502S-E-WDAV

Vials, Screw Top, Glass. Amber, 2ml, with Write-On Patch and fill lines, Deactivated, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Deactivated / Silanized Borosilicate Glass, amberized. Features a wide mouth for easy filling. Low cost for budget minded laboratories. Write-On Patch and fill lines are epoxy paint for reduced smearing with solvents.

Amber vials are often used for light sensitive compounds and considered low actinic.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$41.31

Printed On 2024/06/01