

Catalog Number

9502S-WAV

Vials, Screw Top, Glass. Amber, 2ml, with Write-On Patch and fill lines, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. MicroSolv Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A pharmaceutical grade, Borosilicate Glass, amberized. Features a wide mouth opening for easy filling and ND9 threads. Write-On Patch and fill lines are epoxy paint for reduced smearing.

Amber vials are often used for light sensitive compounds and considered low actinic.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$32.79

Printed On 2024/05/20