

Catalog Number

9502S-WCV

Vials, Screw Top, Glass. Clear, 2ml, with Write-On Patch and fill lines, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. MicroSolv Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A pharmaceutical grade, Borosilicate Glass. Features a wide mouth opening for easy filling and ND9 threads. Write-On Patch and fill lines are epoxy paint for reduced smearing.

This is our Most Popular standard, Borosilicate Glass 2ml Vial. For the Most Popular RSA™ Vial, Click [HERE](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$30.55

Printed On 2024/05/20