

## Catalog Number

9502S-WDCV

**Vials, Screw Top, Glass.** Clear, 2ml, with Write-On Patch and fill lines, Deactivated, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand. 100/PK.



### Additional Info

First hydrolytic, Class A, Type N51A pharmaceutical grade, Borosilicate Deactivated / Silanized Glass. Write-On Patch and fill lines are epoxy paint for reduced smearing. Features a wide mouth opening for easy filling and ND9 threads.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$65.59**

Printed On 2024/05/20