

Catalog Number

9503S-05FB

Inserts, Low Volume for Vials, Glass. Clear, 300ul fill volume with a flat bottom and 6x31mm outer dimensions. For use in 2ml autosampler vials. BASIK Brand. (10 packs of 100). 1000/CS.



Additional Info

First hydrolytic, Class A, Borosilicate Glass. Features a flat bottom. Residual volume is 6ul with most autosamplers.

Use with 9mm Screw Top or 11mm Snap Top 2ml, 12x32mm autosampler vials.

Lowest cost for budget minded laboratories.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$379.06

Printed On 2024/05/19