

## Catalog Number

9503S-09-KR

**Caps, Screw Top, Solid, with a Red-Orange Silicone Rubber / TEF liner, in "knurled", closed Black Caps. For use with 10-425mm thread, 12x32mm autosampler vials. MicroSolv Brand 100/PK.**



### Additional Info

Use with wide mouth autosampler vials. The Silicone Rubber / TEF liner provides a very good seal for longer term storage in autosampler vials.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$50.39**

Printed On 2024/05/20