

## Catalog Number

9503S-35C-K-M

**Caps, Screw Top**, with Silicone Rubber / PTFE Pre-Slit, 40mil Septa, bonded in wide mouth, "knurled", polypropylene, Black Caps. For use with 10-425mm threads, 12x32mm autosampler vials. BASIK Brand.



### Additional Info

Pre-slit for easy penetration by autosampler needles.

This Septa is slightly thinner than most Septa to enhance the ease of needle penetration. 40mil equals 1.016mm.

Click [HERE](#) for smaller pack quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$491.76**

Printed On 2024/05/19