

Catalog Number

9503S-GR-M

Caps, Screw Top, with Silicone Rubber / PTFE Septa fitted in "knurled", polypropylene, Black Caps. For use with 10-425mm thread, 12x32mm wide mouth autosampler vials. Compare to Shimadzu Caps. BASIK



Additional Info

Flange at the bottom of the Cap allows for the vial and cap to be suspended in older Shimadzu autosamplers.

Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$315.56

Printed On 2024/05/19