

Catalog Number

9504R-01

Vial Rack, for 15x45mm Vials. Acrylic, clear, with 40 cavities. 1 EA..



Additional Info

For easy handling and transportation of 4ml, 15x45mm autosampler vials. The clear materials makes filling easier with greater visibility compared to opaque racks. Also, racks are stackable due to the silicone feet on the base. More stable than other racks on the market due to the solid, one piece construction.

Flat bottom vials can be accommodated in this 175.8mm wide x 115.5mm deep and 20mm high rack. Wells are 15.1mm in diameter.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$165.54

Printed On 2024/05/20