

## Catalog Number

9504S-01C-BK

**Caps, Screw Top, Solid, without a liner. "knurled", polypropylene, Black Caps.** For use with 13-425mm thread, 15x45mm autosampler vials.  
MicroSolv Brand 100/PK.



### Additional Info

This cap is not recommended for storage or transportation of chemicals or samples to be analyzed. Good for storage of dry objects or after inserting a liner.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$40.33**

Printed On 2024/05/19