

Catalog Number

9504S-1CP-AQA

Vial & Cap Kit Includes 100 EA 4ml, Screw Top, Clear Autosampler Vials and 100 EA matching Screw Caps with Silicone Rubber / PTFE Septa, bonded in the Black Caps. Vials and Caps are Pre-assembled. MicroSolv



Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps.

Click on links below for more Specifications for the 13mm Thread Vials and Caps:

Supplied with Vials: Cat. No [9504S-0CV](#) and Caps: Cat No. [9504S-10C-BK](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$129.43

Printed On 2024/05/20