

## Catalog Number

9504S-3WAP

**Vial & Cap Kit** Includes 100 EA 4ml, Screw Top, Amber Autosampler Vials with Write-On Patch and fill lines and 100 EA matching Screw Caps with White Silicone Rubber / Tan PTFE Pre-Slit, Septa, bonded in the



### Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps.

Click on links below for more Specifications for the 13mm Thread Vials and Caps:

Supplied with Vials: Cat. No [9504S-WAV](#) and Caps: Cat No. [9504S-30C-BK](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$134.84**

Printed On 2024/06/01