

Catalog Number

9504S-E-WAV

Vials, Screw Top, Glass. Amber, 4ml, with Write-On Patch and graduations, 13-425mm threads, 15x45mm outer dimensions. For use as an autosampler vial. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Borosilicate Glass, amberized. Features ND13 threads. Write-On Patch is epoxy paint for reduced smearing.

Amber vials are often used for light sensitive compounds and considered low actinic.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$42.16

Printed On 2024/05/19