

## Catalog Number

9508C-10C-D

**Caps, Push Top**, with Silicone Rubber / PTFE Septa, fitted in smooth, polyethylene, Natural Transparent Caps. For use with 8mm, Crimp Top autosampler vials. MicroSolv Brand. 100/PK.



### Additional Info

This cap is similar to a Snap Cap but does not need the Snap Ring on a vial to create a tight seal.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$55.89**

Printed On 2024/05/19