

Catalog Number

95096-08-CAP-M

Caps, PE Plug Style, Polyethylene, Natural Transparent Caps. For use with 700ul and 1.0mL, 96 well autosampler plate inserts. AQ Brand. 1000/CS.



Additional Info

Click [HERE](#) for page with compatible inserts. Use 96 Well Plate Insert FILTER.

Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$113.45

Printed On 2024/05/19