


Catalog Number

9509C-30CB

Caps, Snap Top, with Silicone Rubber / PTFE, ultra-pure, Pre-Slit Septa, fitted in smooth, polyethylene, Light Blue Caps. For use with 11mm, 12x32mm, Snap Top autosampler vials. AQR Brand. 100/PK.



Additional Info

This specially designed septa provides easy penetration by autosampler needles with improved resealing compared to other Snap Caps. The ultra-pure Silicone Rubber and PTFE have no added color producing the highest purity septa on the market with virtually no contamination or extractables.

Recommended for LCMS/MS.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$0.00

Printed On 2024/05/20