


Catalog Number

9509S-10C-B

Caps, Screw Top, with Silicone Rubber / PTFE, ultra-pure Septa, fitted in "knurled", polypropylene, Light Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. AQR Brand. 100/PK.



Additional Info

When highest purity in a Septa is required. No color added to the rubber or Teflon for improved purity. Good choice for LCMS, GCMS & LCMS/MS.

For Waters® Alliance, Acquity™ and all other autosamplers that require a Pre-Slit Septa use catalog# [9509S-30C-B](#).

Click [HERE](#) for case quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$0.00

Printed On 2024/05/19