


## Catalog Number

9509S-30B

**Caps, Screw Top**, with White Silicone Rubber / White PTFE, ultra-low bleed, Pre-Slit Septa, bonded in "knurled", polypropylene, Dark Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. AQ



### Additional Info

Pre-slit for easy penetration by autosampler needles and *recommended for Waters® Alliance, Acquity™ and all other autosamplers that require Pre-Slit Septa.*

Good choice for LCMS & LCMS/MS.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$0.00**

Printed On 2024/05/19