

## Catalog Number

9509S-30C-B

**Caps, Screw Top**, with Silicone Rubber / PTFE, ultra-pure, Pre-Slit Septa, fitted in "knurled", polypropylene, Light Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. AQR Brand. 100/PK.



### Additional Info

Highest purity Septa on the market, with no color added to the Rubber or Teflon for improved purity. Good choice for LCMS & LCMS/MS.

Pre-slit for easy penetration by autosampler needles and recommended for Waters® Alliance, Acquity™ and all other autosamplers that require Pre-Slit Septa. For lower cost, non-slit version please click [HERE](#) for more information.

Click [HERE](#) for case quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$0.00**

Printed On 2024/05/19