

## Catalog Number

9509S-3WCP-RS

**Vial & Cap Kit** Includes 100 EA 2ml, Screw Top, Clear RSA™ Autosampler Vials with Write-On Patch and fill lines and 100 EA matching Screw Caps with Pre-Slit, Clear AQR Silicone Rubber / Clear



### Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use. This is our Most Popular Easy Purchase Pack for RSA™ Vials with Pre-Slit Septa.

Click on links below for more Specifications for the 9mm Vials and Caps:

Supplied with Reduced Surface Activity Vials: Cat. No [9509S-WCV-RS](#) and LCMS/MS compatible Caps: Cat No. [9509S-30C-B](#).

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$0.00**

Printed On 2024/05/19