

## Catalog Number

9512C-0-0AV-P

**Vials, Crimp Top, Glass. Amber, 1.5ml, U Shape. An 11mm Serum Finish with 12x32mm outer dimensions and a fluted stand. For use as an autosampler vial. MicroSolv Brand. 100/PK.**



### Additional Info

First hydrolytic, Class A, Type N51A, Borosilicate Glass, amberized. Amber vials are often used for light sensitive compounds and considered low actinic.

Compares to Champagne Style Vials. Features a 15ul residual volume for most autosamplers due to the center draining, inner V shape.

*Use with 11mm aluminum crimp caps. These lower cost center draining vials have serum finish tops that provide the best seal for chromatography, when used properly with correct aluminum crimp caps.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$79.06**

Printed On 2024/05/30