

## Catalog Number

9512C-0-0DCV

**Vials, Snap Top, Glass. Clear, 1.8ml, MRQ Style. Deactivated. An 11mm snap-ring and 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand. 100/PK.**



### Additional Info

First hydrolytic, Class A, Type-N51A, Deactivated / Silanized, Borosilicate Glass. The MRQ Style Vials have a fill volume of 1.2ml and a very low residual volume of 5ul due to the center-draining effects and the conical bottom with a well that is slightly larger than autosampler vial needles. Side entry port needles will recover the maximum amount during use.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$267.85**

Printed On 2024/05/19