

Catalog Number

9512C-0AV-P

Vials, Snap Top, Glass. Amber, 1.5ml, U Shape. An 11mm snap-ring with 12x32mm outer dimensions and a fluted stand. For use as an autosampler vial. MicroSolv Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A, Borosilicate Glass.

Compares to Champagne Style Vials. Features a 15ul residual volume for most autosamplers due to the center draining, inner U shape.

We recommend to use these vials with [Soft-Guard™](#) Snap Caps.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$83.16

Printed On 2024/05/20