

Catalog Number

9512C-0E-0CV

Vials, Crimp Top, Glass. Clear, 1.8ml, Max Recovery Style. An 11mm Serum Finish with 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. MicroSolv Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A, Borosilicate Glass.

Features an 25ul residual volume for most autosamplers due to the center draining interior and conical bottom.

Use with 11mm aluminum crimp caps. These lower cost center draining vials have serum finish tops that provide the best seal for chromatography, when used properly with correct aluminum crimp caps.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$89.04

Printed On 2024/05/19