

Catalog Number

9512C-1MP-RS

Vial & Cap Kit includes 100 EA 1.2ml, MRQ, Snap Top, Clear RSA™ Autosampler Vials and 100 EA matching Snap Caps with White Silicone Rubber / Tan PTFE Septa, fitted in the Blue Caps. RSA Brand Easy



Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the 11mm Vials and Caps:

Supplied with Reduced Surface Activity Vials: Cat. No [9512C-0CV-T-RS](#) and Caps: Cat No. [9502C-11-MB](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price
\$0.00

Printed On 2024/06/01