

Catalog Number

9512C-AQR-1CP

Vial & Cap Kit includes 100 EA 1.2ml, MRQ, Snap / Crimp Top, Clear Autosampler Vials and 100 EA matching Crimp Caps with Clear AQR Silicone Rubber / Clear PTFE, ultra-pure Septa, fitted in the Silver Caps.



Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the Crimp Top Vials and Caps:

Supplied with Vials: Cat. No [9512C-0CV-T](#) and Caps: Cat No. [9509C-AQR-1S](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$223.34

Printed On 2024/05/20