

Catalog Number

9512S-1AP-RS

Vial & Cap Kit Includes 100 EA 1.8ml, Maximum Recovery, Screw Top, Amber RSA™ Autosampler Vials and 100 EA matching Screw Caps with Clear AQR Silicone Rubber / Clear PTFE, ultra-pure Septa, fitted in the



Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the 9mm Vials and Caps:

Supplied with Reduced Surface Activity Vials: Cat. No [9512S-0AV-RS](#) and LCMS/MS compatible Caps: Cat No. [9509S-10C-B](#).

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$275.26

Printed On 2024/05/20