

## Catalog Number

9512S-KDAP-PX

**Vial & Cap Kit.** 100 EA 1.8ml, Screw Top, Hydrophobic, Amber Autosampler Vials and matching Solid Black Caps with Clear Silicone Rubber / PTFE Liner. RSA-Pro X Brand. 1 PK.



### Additional Info

These 1.8ml RSA, Glass, Max Recovery Autosampler Vial & Cap Kit Packs (*Easy Purchase Packs*) are supplied with [Amber Center Draining RSA Vials](#) that have been treated with [RSA-Pro X™ Technology](#) to prevent their already pristine surface from adsorbing stubbornly sticky compounds that have affinity for the glass structure. The Hydrophobic surface is not only non-adsorbing for most compounds but is extremely stable even in the presence of Water. The hydrolytic stability will last for years compared to Silanized Vials that will usually degrade in the presence of Water. This product is manufactured in an ISO 9001:2015 Certified facility.

The supplied, Solid Top matching [Black, AQR™ Screw Caps](#) have the most inert Silicone Rubber / PTFE Liner available; for total Vial and Cap inertness. Often chosen for storage or shipping of standards, proteins, peptides, enzymes, natural products and other precious compounds.

*These treated vials are only sold in RSA-Pro X™ Kits. Click [HERE](#) for replacement Caps.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$482.16**

Printed On 2024/05/19