

## Catalog Number

95296-LG-30R-EA

**Mat, Sealing, Large, Chemical Resistant Polymer, Pre-Slit, Plug Style for 96, 500 or 1000ul Inserts. U-2D Brand. 1 EA.**



### Additional Info

Mats are made of a chemically resistant polymer that can be used with typical HPLC organic solvents. Designed and molded to seal up to 96 each 500ul or 1ml Inserts in the U-2D Micro-Sample Management System. These mats are Pre-Slit and have alphanumeric markings on the top side of the mats, along the edges for identification of the Inserts.

Click [HERE](#) for case quantity.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$26.09**

Printed On 2024/05/20