

## Catalog Number

9532C-0-WAV

**Vials, Snap Top, Glass. Amber / Clear Insert, 300ul with Write-On Patch, Fused Insert Style. An 11mm Serum Finish with 12x32mm outer dimensions. For use as an autosampler vial. AQ Brand. 100/PK.**



### Additional Info

Vial and inserts are made with first hydrolytic, Class A, Type N51A, Borosilicate Glass, amberized. The Write-On Patch is epoxy paint and will resist smearing when working with solvents. This product is manufactured in an ISO 9001:2015 Certified facility.

*Amber vials are often used for light sensitive compounds and considered low actinic.*

*Use with 11mm aluminum crimp caps. These vials have serum finish tops that provide the best seal for chromatography, when used properly with correct aluminum crimp caps.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$204.20**

Printed On 2024/05/19