

## Catalog Number

9532C-0-WCV

**Vials, Crimp Top, Glass. Clear / Clear Insert, 300ul, with Write-On Patch and fill lines, Fused Insert Style. An 11mm Serum Finish with 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand.**



### Additional Info

Vials and Inserts are made with first hydrolytic, Class A, Type-N51A Clear, Borosilicate Glass. The insert has a 300ul total fill volume and a 3ul residual volume for many autosamplers and autosampler needles. The Write-On Patch and fill lines are epoxy paint and will resist smearing when working with solvents.

*Use with 11mm aluminum crimp caps. These lower cost Fused Insert vials have serum finish tops that provide the best seal for chromatography, when used properly with correct aluminum crimp caps.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$301.06**

Printed On 2024/05/19