

Catalog Number

9532C-TS-C

Vials, Snap Top, Plastic. Clear / Clear Glass Insert, 200ul Fused Insert™. An 11mm snap-ring and 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand. 100/PK.



Additional Info

These narrow mouth vials are made with TPX (polymethylpentene) Plastic and have a first hydrolytic, Type N51A, Clear, pharmaceutical grade, Borosilicate Glass Insert, with a conical bottom, embedded in the plastic. These vials have a 200ul total fill volume and will not leak from the bottom if the needle crashes the insert. The Insert will not move in the vial.

Compares to TopSert, GlassPlastic™ and Glastic™ Brands.

Click [HERE](#) for Amber Vials.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$246.74

Printed On 2024/05/19