

Catalog Number

9534S-1DCP

Vial & Cap Kit Includes 100 EA 3.8ml, Maximum Recovery, Screw Top, Clear Deactivated Library Storage Vials and 100 EA matching Screw Caps with Silicone Rubber / PTFE Septa bonded in the Black Caps. AQ



Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps.

Click on links below for more Specifications for the 13mm Vials and Caps:

Supplied with Silanized / Deactivated Vials: Cat. No [9534S-0DCV](#) and Caps: Cat No. [9504S-10C-BK](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$377.65

Printed On 2024/05/20