


## Catalog Number

9552S-10C-W

**Caps, Screw Top**, with Silicone Rubber / PTFE Septa, bonded in "knurled", polypropylene, White Caps. For use with 16mm thread, 22x37mm autosampler vials. MicroSolv Brand. 100/PK.



### Additional Info

Recommended for use with the 5mL Center Draining vials. Click [HERE](#) for compatible vials.

*Compared to the caps used in the Agilent InfinityLab™ LC Series; fully compatible.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$96.74**

Printed On 2024/05/19