

Catalog Number

9572S-PP-K

Vials, Screw Top, Plastic. 700ul, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand. 1000/CS.



Additional Info

These injection molded, translucent, plastic vials are made with virgin polypropylene and chosen for bio-analytical reasons.

Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$302.60

Printed On 2024/05/20