

Catalog Number

9582S-0AV

Vials, Screw Top, Glass. Amber, 2ml, 8-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Borosilicate Glass, amberized. Features a narrow mouth which may be useful to reduce evaporation and requires an 8mm screw cap.

Amber vials are often used for light sensitive compounds and considered low actinic.

Lowest cost in class.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$33.76

Printed On 2024/05/20