

## Catalog Number

9582S-0DCV

**Vials, Screw Top, Glass.** Clear, 2ml, Deactivated, 8-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. BASIK Brand. 100/PK.



### Additional Info

First hydrolytic, Class A, Deactivated / Silanized, Borosilicate Glass. Features a narrow mouth which may be useful to reduce evaporation and requires an 8mm screw cap.

Good choice for biological compounds such as proteins and peptides due to the hydrophobic nature of the vial.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$60.56**

Printed On 2024/05/20