

## Catalog Number

9582S-30C-K

**Caps, Screw Top**, with 1.3mm thick Silicone Rubber / PTFE Pre-Slit Septa, fitted in "knurled", polypropylene, Black Caps with 5.5mm center holes. For use with 8-425mm thread, 12x32mm autosampler vials.



### Additional Info

Pre-Slit for easy penetration by autosampler needles.

Click [HERE](#) for case quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$69.24**

Printed On 2024/05/19