

## Catalog Number

9582S-WCV

**Vials, Screw Top, Glass.** Clear, 2ml, with Write-On Patch and fill lines, 8-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. BASIK Brand. 100/PK.



### Additional Info

First hydrolytic, Class A, Borosilicate Glass. Write-On Patch and fill lines are epoxy paint which will resist smearing during use. Features a narrow mouth which may be useful to reduce evaporation.

The narrow opening and write-on patch make these ideal for stability testing and sample storage.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$33.90**

Printed On 2024/05/20