

## Catalog Number

97010-CV05

**Vials, Push Top, Storage, Glass, Clear, 5ml. Economical, with ND18 tops. For use with Econo-Storage Push Top Caps. BASIK Brand. 100/PK.**



### Additional Info

These Vials are made with 3rd Hydrolytic Glass and quickly open and reseal due to the "Snap" design. These 5ml Vials have a total fill volume of 9.1ml with 40x20mm outer dimensions. Not made for Crimp Caps.

Use [97010-SC18](#) Caps with these Vials, sold separately. These vials are not intended for Analytical Testing and should only be used for Storage.

*Low cost for budget minded laboratories.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$32.91**

Printed On 2024/05/19