

Catalog Number

97018-CV-16

Vials, Screw Top, Storage, Glass, Clear, 16ml, with 18-400mm Screw Threads. For use with For use as a Storage Vial. MicroSolv Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A, Pharmaceutical Grade, Borosilicate Glass.

These Storage Vials have 71x20.6mm outer dimensions and a total fill volume of 17.4ml and can be used for long term storage with a proper cap.

Use with Caps that have 18-400 or ND18 screw threads, sold separately.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$89.18

Printed On 2024/05/19