

Catalog Number

97031-0CV

Vials, EPA Style, Glass. Clear, 30ml. 70x27.5mm outer dimensions. For use as an EPA test vial or storage bottle. ND24 threads. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type 33, Borosilicate Glass. Fill volume is approximately 31ml.

Use 24mm screw caps that have 24-414 threads, sold separately. *Lower cost for budget minded laboratories.*

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$79.31

Printed On 2024/05/31