

## Catalog Number

97040-05C-K

**Caps, Screw Top**, with White Silicone Rubber / White PTFE Septa, fitted in knurled polypropylene, Black Caps. For use with 24-414mm threaded EPA Style vials. BASIK Brand. 100/PK..



### Additional Info

EPA quality septa can be used for EPA testing or for long term storage with any 30, 40 or 60 EPA Style vial.

Cap has a 12.5mm top hole.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$67.01**

Printed On 2024/05/20