

Catalog Number

97040-0AV

Vials, EPA Style, Glass. Amber, 40ml. 95x27.5mm outer dimensions. For use as an EPA test vial or storage bottle. ND24 threads. MicroSolv Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A, Pharmaceutical grade Borosilicate Glass, amberized. Amber vials are often used for light sensitive compounds and considered low actinic. Fill volume is approximately 42.9ml.

Use 24mm screw caps that have 24-400 threads, sold separately.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$200.96

Printed On 2024/05/19