

## Catalog Number

PQA-L06-250

**Linearity L6 Solution.** For UV-vis qualification. Chemical Solutions brand. 250 ml. 1 EA.



### Additional Info

Use to test your HPLC or UV-vis detector for linearity precision and accuracy with the Chemical Solutions PQ or HSQ kits. This replacement solution is nominally labelled as Linearity L6 Solution and is NIST traceable. Supplied in a 250 mL plastic bottle and leak proof cap.

Click [HERE](#) for specifications and concentration of this solution.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$485.21**

Printed On 2024/05/19