

## Catalog Number

04071-C25

**Capillary, Bare Fused Silica, 75um ID by 25 meters. Simplus Brand. 1 EA.**



## Additional Info

This untreated, fused silica capillary is clad in polyimide and designed for use in Capillary Electrophoresis. Supplied on a plastic, non static spool for storage.

Manufactured with extreme precision, these capillaries have excellent concentricity throughout the entire length of the capillary. The outer diameter (OD) is 365um. Made in USA.

Click [HERE](#) for 1 meter length and [HERE](#) for 10 meter length.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$821.54**

Printed On 2024/05/19