

## Catalog Number

07220-S

**Syringe, Replacement,** for Hand Held Low Volume Pump, 500ul Volume, 365um OD Tubing. MicroSolv Brand. 1 EA.



### Additional Info

This plastic Syringe is a Replacement accessory for the MicroSolv Hand Held, Low Volume Pump for 365um OD CE or CEC Capillaries. Supplied with 365um ID low pressure Plastic Tubing, that slips over the Syringe Needle and Capillary with a negative tolerance.

Click [HERE](#) for Hand Held Pump and [HERE](#) for Replacement Tubing ordering information.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$15.73**

Printed On 2024/05/20