

## Catalog Number

40000-05P-2

**HPLC Column, Silica-C**, 4um, 2.1mm ID x 50mm Length, 100A. Cogent  
TYPE-C Brand. 1 EA.



### Additional Info

The Silica-C™ HPLC Columns are made with [Cogent TYPE-C™ Silica](#), and provide unique Selectivity due to its surface chemistry (especially in Normal Phase HPLC) compared to standard silica, This Silica Hydride Phase is not modified and therefore has no end-capping or organic ligands.

Chosen for analyzing compounds that are not soluble in Water including but not limited to: lipids, peptides, proteins, oligonucleotides and other compounds that require non polar solvents.

- Faster Equilibration than other Silica Columns
- No need to control Moisture in Mobile Phase Solvents
- Robust Column Lifetimes

*This short and narrow-bore column provides short run times but requires a pump that can precisely deliver, and a flow cell that made for, reduced flow rates with lower volumes. This*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$839.68**

Printed On 2024/05/19