

## Catalog Number

40000-P50-50

**Prep HPLC Column, Silica-C, 10um, 100A. 50mm ID x 50mm Length, Preparative. Cogent TYPE-C Brand. 1 EA.**  
**CUSTOM ORDER.**



### Additional Info

This Cogent Silica-C™ Preparative, HPLC Column is made with 4um, 100A, [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent yields and turn around times.

Chosen for analyzing compounds that are not soluble in Water including but not limited to: lipids, peptides, proteins, oligonucleotides and other compounds that require non polar solvents.

- Faster Equilibration than other Silica Columns
- No need to control Moisture in Mobile Phase Solvents
- Robust Column Lifetimes
- Use Minimal Salt to Elute in ANP.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$22,954.86**

Printed On 2024/05/19