

## Catalog Number

40031-P21-50

**Prep HPLC Column, UDA (wcx), 4um, 100A. 21.2mm ID x 50mm Length, Preparative. Cogent TYPE-C Brand. 1 EA.**



### Additional Info

This Cogent UDA™ (wcx) HPLC Column is made with [Cogent TYPE-C™ Silica](#), and has an 11 carbon chain with a terminal carboxylic acid group. This silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Extremely Robust Column Lifetimes.

Use in Reversed Phase (RP), Aqueous Normal Phase (ANP) or as a Weak Cation Exchange Column (wcx at pH 6 or higher).

- Faster Equilibration than other Silica Columns
- No need to control Moisture in Mobile Phase Solvents
- Robust Column Lifetimes
- Use Minimal Salt to Elute in ANP.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$4,884.58**

Printed On 2024/05/18